




T1BBDCD456U2



A2001 A2001 A2401



A2601 A2501 A2621



Blood Bags

- Improved plastic is ideal for storage of platelets
- XT-612's superb gas permeability allows platelets to be stored for 5 days with good viability
- The Terumo lancet point needle minimises donor discomfort
- Packed flat in Terumo's exclusive blister package - lay flat for quick and easy removal
- Each style of blood bag is packed in a specially designed, colour coded, blister package for convenient storage



Product Code	Description	Unit
T1BBDCD456U2	Double - Additive, 450ml	10/Box
T1BBSCD456U	Single - Additive, 450ml	5/Box
T1BBT060CBBL	Transfer, 600ml	Each
T1BBT060CBB	Transfer, 600ml	25/Box
T1BBT100BB71	Transfer, 1000ml	10/Box
T1BBT200BB71	Transfer, 2000ml	25/Box

Blood Collection Tubes

Flat Bottom, Plastic



Product Code	Description	Width x Height	Unit
A2001	2.5ml With K2 Edta - Green	16 x 60mm	50/Box
A2201	2.5ml With Kf+Na2 Edta - Orange	16 x 50mm	50/Box
A2401	2.5ml With Li Heparin Edta - Blue	16 x 60mm	50/Box
A2407	10ml With Li Heparin Edta - Dark Blue	16 x 100mm	50/Box
A2601	0.25ml With Na Citrate Pink Cap - Pink/Yellow	16 x 50mm	50/Box
A2501	0.4ml With Na Citrate Yellow Cap - Yellow	16 x 50mm	50/Box
A2621	0.5ml With Na Citrate Pink Cap - Pink	16 x 50mm	50/Box

Venoject II Evacuated Plastic

- Safety features: Shatter-resistant, plastic evacuated blood collection tubes lessen risk of exposure to blood
- UltraSeal™ closure reduces potential for splatter and biohazard contamination
- Efficiency in clinical processing: Clear plastic evacuated blood collection tubes virtually eliminate breakage in centrifuge
- UltraSeal™ tube closure minimises red cell hang-up
- In chemistry, time-saving clot activator disc promotes quicker clotting
- Piercing sampler allows for slide preparation in seconds without opening tube
- Venoject® II Tubes incinerate cleanly with no poisonous residues
- In haematology and coagulation, smooth plastic surface helps prevent clotting



Product Code	Description	Capacity	Unit
TVP053SDK	Purple - Edta K2	2.5ml	100/Box
TVP050SFX	Grey - Fluoride/Oxalate	5ml	100/Box
TVF054SHLW	Green - Lithium Heparin	5ml	100/Box
TVF054SPW	Plain Silicone Coated	5ml	100/Box
TVP076SAS	Cherry Red - Autosep Serum Gel	7ml	140/Box
TVP100SP	Red - Plain Silicone Coated	10ml	100/Box

Blood Lancets

Medipoint



Product Code	Description	Unit
JPBLCT03	Point Length 3.7mm	200/Box
JPBLCT04	Point Length 3.7mm with Plastic Encapsulated	200/Box

- GST Exempt Item



Product Code	Description	Unit
BA3030	Minilet Blood Lancets Short Point	200/Box

- GST Exempt Item



Product Code	Description	Unit
UNISTIK	Unistik Blood Letting Devices	50/Box

- GST Exempt Item



