

## Examination Couch

- Adjustable back rest with comfortable seamless vinyl upholstered top
- Finish: Epoxy powder coated with vinyl cover
- Safe working load: 130kg



Product Code	Description	Unit
KEKH320G	1,800 x 600 x 850mm	Each

- Extra Freight Charges Apply



## Examination Lighting

### HL 1200 Light

- Built-in temperature filter - Avoids patient discomfort
- Swivels through complete 180°. Flexible and easy to position.
- 'Stay-Put' construction of flexible section - Adjust the light to virtually any position and it 'stays put', will not sag or drop
- 20 W IRC reflector lamp with innovative bulb technology
- Long Life 12V/20W IRC reflector lamp, 4000 hrs
- More than twice the brightness in comparison with conventional 20 Watt lamps for brightness
- 90 000 lux light output at 30cm. Bright, white light



### Wall Mounted

Product Code	Description	Unit
H0527052	With wall-mounting bracket can also be fitted to optional wheeled stand	Each

### With Mobile base

Product Code	Description	Unit
H0527068	With wheeled stand, metal base	Each
H0527078	With lightweight stand, Polycarbonate base	Each
H0527056	With additional multi-purpose clamp	Each
H0527075	Spare bulb 12V / 20W IRC reflector lamp	Each

- Extra Freight Charges Apply



## Operating/Magnification Lighting

- Shadow free light and quality magnification
- Assembled with 4 Diopter Lens, providing magnification of double actual size



Product Code	Description	Unit
LUXOMAGLTW	LUXO K - Light & Magnifier	Each
LUXOWALLBKT	Wall Bracket B	Each
LUXOSTAND	Mobile Stand	Each
LUXOTUBE	Fluorescent Light Tube	Each

- Extra Freight Charges Apply



## Foot Stools & Steps



Product Code	Description	L x H x W	Unit
KEKH32SP	Single Step	500 x 200 x 250mm	Each
KEKH32DP	Double Steps	500 x 300 x 250mm	Each

- Extra Freight Charges Apply



## Privacy Screen

### 3 Bay

- Easy folding, mobile, light weight, heavy duty screen with plastic hinge and low ground clearance for privacy.
- Finish: Epoxy powder coated



Product Code	Description	Unit
KEKH514G	3450 x 550 x 1800 mm (L x W x H)	Each

- Extra Freight Charges Apply





## Trolley

### Dressing

- Stainless steel shelves and epoxy powder coated side frames
- Comes with push handle on both sides and guard rail on 3 sides



Product Code	Description	L x H x W	Unit
KEKH140	No Drawer	500 x 500 x 1,015mm	Each
KEKH140D	With Drawer	500 x 500 x 1,015mm	Each
KEKH150D	With Drawer	800 x 500 x 1,050mm	Each

- Extra Freight Charges Apply



**Livingstone**  
**offers a complete range of Sutures**  
 For our complete product range please  
 contact us on **1300 780 078**.