

• Extra freight charges will apply for all hazardous material handling equipment



Storage Cabinets, Corrosive



Blue Powder Coated, Metal with Polyethylene Trays

• Meets Australian standard AS 3780

Compact



Product Code	Description	Dimension (H x W x D)	Unit
SSC25714B	30L, 2 Shelves, 1 Door	800 x 595 x 465mm	Each
SSC25712B	60L, 2 Shelves, 1 Door	1,005 x 595 x 465mm	Each



Under Bench

Product Code	Description	Dimension (H x W x D)	Unit
SSC25748B	100L, 2 Shelves, 2 Door	800 x 920 x 565mm	Each



Standard

Product Code	Description	Dimension (H x W x D)	Unit
SSC25302B	160L, 2 Shelves, 2 Door, 6 x 20L Drums	1,230 x 1,100 x 465mm	Each
SSC25452B	250L, 3 Shelves, 2 Door, 9 x 20L Drums	1,760 x 1,100 x 465mm	Each



Blue High Density Polyethylene



Product Code	Description	Dimension (H x W x D)	Unit
SSC24040	15L, 1 Door	500 x 365 x 415mm	Each
SSC24160	30L, 2 Shelves, 2 Door	890 x 920 x 565mm	Each

Storage Cabinets, Flammable



Yellow Powder Coated, Metal

• Meets Australian standard AS 1940



Compact

Product Code	Description	Dimension (H x W x D)	Unit
SSC25714	30L, 2 Shelves, 1 Door	800 x 595 x 465mm	Each
SSC25712	60L, 2 Shelves, 1 Door	1,005 x 595 x 465mm	Each



Under Bench

Product Code	Description	Dimension (H x W x D)	Unit
SSC25748	100L, 2 Shelves, 2 Door	800 x 920 x 565mm	Each



Standard

Product Code	Description	Dimension (H x W x D)	Unit
SSC25302	160L, 2 Shelves, 2 Door, 6 x 20L Drums	1,230 x 1,100 x 465mm	Each
SSC25452	250L, 3 Shelves, 2 Door, 9 x 20L Drums	1,760 x 1,100 x 465mm	Each
SSC25602	250L, 3 Shelves, 2 Door, 12 x 20L Drums	1,760 x 870 x 870mm	Each



Drum

Product Code	Description	Dimension (H x W x D)	Unit
SSC25662	207L, Horizontal Drum, 2 Door	1,385 x 765 x 1,225mm	Each
SSC25701	207L, Vertical Drum, 2 Door with Rollers	1,760 x 870 x 870mm	Each

Storage Cabinets, Organic Peroxide

Yellow Powder Coated, Metal

- Meets Australian standard AS 2714

Compact



Product Code	Description	Dimension (H x W x D)	Unit
SSC25712OP	50L, 2 Shelves, 1 Door	1,005 x 595 x 465mm	Each

Under Bench

Product Code	Description	Dimension (H x W x D)	Unit
SSC25748OP	60L, 2 Shelves, 2 Door	800 x 920 x 565mm	Each

Standard

Product Code	Description	Dimension (H x W x D)	Unit
SSC25452OP	100L, 3 Shelves, 2 Door, 9 x 20L Drums	1,760 x 1,100 x 465mm	Each



Under Bench

Product Code	Description	Dimension (H x W x D)	Unit
SSC25748T	100L, 2 Shelves, 2 Door	800 x 920 x 565mm	Each

Standard

Product Code	Description	Dimension (H x W x D)	Unit
SSC25302T	160L, 2 Shelves, 2 Door, 6 x 20L Drums	1,230 x 1,100 x 465mm	Each
SSC25452T	250L, 3 Shelves, 2 Door, 9 x 20L Drums	1,760 x 1,100 x 465mm	Each

Storage Cabinets, Oxidising Agents

Yellow Powder Coated, Metal

Meets Australian standard AS 4326

Compact



Product Code	Description	Dimension (H x W x D)	Unit
SSC25712OA	60L, 2 Shelves, 1 Door	1,005 x 595 x 465mm	Each

Under Bench

Product Code	Description	Dimension (H x W x D)	Unit
SSC25748OA	100L, 2 Shelves, 2 Door	800 x 920 x 565mm	Each

Standard

Product Code	Description	Dimension (H x W x D)	Unit
SSC25452OA	250L, 3 Shelves, 2 Door, 9 x 20L Drums	1,760 x 1,100 x 465mm	Each

Shelves for Cabinets

Metal



Product Code	Description	Unit
SSC25901	Shelf to Suit SSC25302 and SSC25452 Cabinets	Each
SSC25909	Shelf to Suit SSC25602 Cabinets	Each
SSC25995	Shelf to Suit SSC25714 and SSC25712 Cabinets	Each

Storage Cabinets, Pesticide

Green Powder Coated, Metal

- Meets Australian standard AS 4452
- A safety shower and eye wash must be installed where packages are opened



Under Bench




Product Code	Description	Dimension (H x W x D)	Unit
SSC25748P	100L, 2 Shelves, 2 Door, 10 x 10L Drums	800 x 595 x 920mm	Each

Standard

Product Code	Description	Dimension (H x W x D)	Unit
SSC25302P	160L, 2 Shelves, 2 Door, 6 x 20L Drums	1,230 x 1,100 x 465mm	Each
SSC25602P	250L, 3 Shelves, 2 Door, 12 x 20L Drums	1,760 x 870 x 870mm	Each
SSC25452P	250L, 3 Shelves, 2 Door, 9 x 20L Drums	1,760 x 1,100 x 465mm	Each

Safety Cans, Dispensing

Red, Terne Plate with 19mm Flexible Hose



Product Code	Description	Unit
SSC10327	3.8L	Each
SSC10527	7.6L	Each
SSC10728	11L	Each
SSC10828	19L	Each

Red, Terne Plate with Bottom Faucet



Product Code	Description	Unit
SSC10307	3.8L	Each
SSC10707	11L	Each
SSC10807	19L	Each

Safety Cans, Disposal

Red, High Density Polyethylene with Plated Steel Fittings



Product Code	Description	Unit
SSC14722	7.6L	Each
SSC14755	19L	Each

- Also available with stainless steel fittings



Storage Cabinets, Toxic Substances (Poison)

White Powder Coated, Metal

Meets Australian standard AS 4452

Compact



Product Code	Description	Dimension (H x W x D)	Unit
SSC25714T	30L, 2 Shelves, 1 Door	800 x 595 x 465mm	Each
SSC25712T	60L, 2 Shelves, 1 Door	1,005 x 595 x 465mm	Each

Safety Cans, Storage

• These safety storage cans and dispensing cans have a built-in flash arrester. The flash arrester stops naked flame from entering the can during either filling or dispensing thereby eliminating any risk from flashback causing explosion. The contents of the can are fully protected from fire should one occur as the flash arrester allows for increased pressure to vent off safely without causing an explosion.



Red, Terne Plate



Product Code	Description	Unit
SSC10001	0.5L	Each
SSC10101	0.9L	Each
SSC10201	1.9L	Each
SSC10301	3.8L	Each
SSC10501	7.6L	Each
SSC10701	11L	Each
SSC10801	19L	Each
SSC11200	Pouring Funnel, Polyethylene	Each



Red, High Density Polyethylene with Plated Steel Fittings



Product Code	Description	Unit
SSC14045	1.9L, Oval	Each
SSC14140	3.8L, Oval	Each
SSC14251	10L, Round	Each
SSC14501	19L, Round	Each
SSC11201	Pouring Funnel, Polyethylene	Each

Red, High Density Polyethylene with Stainless Steel Trim



Product Code	Description	Unit
SSC14561	19L	Each

HPLC Waste Collection Solutions



High Density Polyethylene with Quick Disconnect & Stainless Steel Fittings (Ideal for HPLC Applications)

- Quick disconnect allows the free flow of solvent into container through tubing
- Vapour and fire protection
- Safeguard against accidents
- Maximise safety and convenience
- Manifold systems for multiple units



Product Code	Description	Unit
SSC12752	7.6L	Each
SSC12755	19L	Each

Spill Control Pallets

- Seamless construction - eliminate leaks
- Chemical resistant - polyethylene construction
- Multiple structural supports to handle heavy loads
- Self locating - can be removed easily for clean up
- Wide forklift pockets for easy relocation
- Load Bearing Capacity - 1,134kg (2 drum); 2,268kg (4 drum)



Product Code	Description	Dimension (H x W x D)	Unit
SSC28234	2 Drum Inline	390 x 1,250 x 640mm	Each
SSC28254	4 Drum Square	260 x 1,250 x 1,250mm	Each

Spill Containment Caddy Karts

- Safe drum handling with mobile dispensing capability and in-built containment
- Chemical resistant - polyethylene construction
- 250L capacity with open view sump for easy spill detection and clean up
- Double loop straps hold drum securely in place
- One piece body eliminates leaks



Product Code	Description	Unit
SSC28908	Outdoor (with Pneumatic Tyre)	Each
SSC28901	Indoor (Solid Plastic Wheels)	Each

Safe Cigarette Butt Disposal

- Self extinguishing the secret, weatherproof
- Safe and tidy, insurance approved for safety
- Low maintenance, internal galvanized collection bucket
- Stable, lockable tie down or weighted security



Product Code	Description	Unit
SSC26800B	Beige	Each
SSC26800DB	Black	Each
SSC26800CB	Chocolate Brown	Each
SSC26800FG	Forrest Green	Each
SSC2600PG	Grey	Each
SSC26800Y	High Visibility Yellow	Each
SSC26800T	Terricotta	Each